



360 Ultra HIFU

High Intensity
Focused Ultrasound

360 Ultra HIFU


360 Ultra HIFU, a non-invasive treatment, goes deeper to stimulate production of your body's own collagen and elastin, for natural-looking results without surgery or downtime.

3

Session

||



10 years younger




Theory

What is 360 Ultra HIFU?

High Intensity Focused Ultrasound (HIFU) delivers the heat energy to SMAS layers and subcutaneous tissue that can stimulate and renew the skin and therefore improving the skin texture and reducing sagging skin. This technique can be applied to the face as well as the whole body, and also, it works equally well on every skin tone.



Before the treatment, collagen in our skin is weak and disorganized which makes our skin saggy with fine lines or wrinkles.



Micro focused ultrasonic rays or energy passes through the skin tissue without damaging the skin cells, but stimulating the collagen in the skin.



The skin responds to the stimulations, making up new collagen and elastin at the same time, so the skin becomes tighter and recovers youthful state after one month.

Handpiece


➤ Classic HIFU handle: Fixed-spot operation

As the traditional ultrasonic anti-aging treatment, the pain is strong, while the golden ultrasonic handle continues the linear output mode of the classic ultrasonic technology, and applies the ultra-pulse technology, which has obvious improvement in comfort and energy output speed.

➤ Ultra gun handle: In motion operation



Focusing on a single point, aiming at the treatment needs of dead corners, the ultra gun handle is applicable for in-motion operation and can offer faster treatment. The design of the ultra gun handle complements the difficulty that the classic ultrasonic handle cannot treat dead corners. It eliminates the forbidden zone that the fixed-point operation has, and it is painless, which brings a new breakthrough for ultrasonic anti-aging treatment.




Cartridge

➤ 8 Cartridge

8 cartridges can penetrate different depths of the face and body precisely, and are multifunctional. Small investments yield big returns.

For Facial Treatment



DS10-1.5
(Optional)

DS7-3.0
(Standard)

DS4-4.5
(Standard)



For Facial Treatment



DS10-1.5
(Optional)

DS7-3.0
(Standard)

DS4-4.5
(Standard)



DS2-13.0

DS3-7.0

Standard




Interface

➤ Friendly Software

The new interface is very user-friendly, easy-to-use and succinct.



3 handles systems, clear to choose



professional system, easy to set
parameter



360 Ultra HIFU system, showing exact
treatment areas

Application

➤ Suitable for what kind of conditions ?


Ideal for removing sagging skin, aging phenomenon for any areas such as forehead, eye brows, cheekbones, risorius, cheeks, chin, neck and body.



- Face lifting & Tightening
- Rebuild Risorius
- Improve Skin Tones
- Improve Double Chin
- Wrinkle Removal
- Remove Nasolabial Folds
- Body Slimming & Tightening

Advantage

- **Ultra-pulse Technology**
High clinical treatment acceptance, the skin getting warm as sliding heating technology kicks in, heat accumulation under the skin, no application of anesthesia, good experience, high comfort level.
- **No restricted area**
360° full coverage, more efficient.
- **Layered Anti-aging**
comprehensively covered.
- **SMAS contraction**
Collagen remodeling, elastine fibril contraction.
- **Instant result**
After treatment, you will see 20-30% of the immediate result, and the best result will manifest in 3-6 months.
- **No downtime & Totally non-invasive**
The treated area will become red in the first several hours, and then recover.
- **Year-round**
refined treatment Combination of HIFU handle and ultra gun, achieve full coverage of the treatment area.
- **Lasting effect**
Good result will last 1-2 years at least.



Benefits



HIFU can generate up to
\$500/session

Treatment areas	Charge reference	Sessions	Revenue
Face	\$300	2	\$600
Arm	\$200	2	\$400
Abdomen	\$300	2	\$600
Thigh	\$300	2	\$600
Calf	\$200	2	\$400

Testimonials



Dr Justin Boey Medical Director
Beauty Clinic Owner

My clients like the treatment very much. They enjoy the process very much. I am very happy to see my clients skin get better.



Dr Haus Ranu

Our clinic uses 360 Ultra HIFU machine for a while, it has the same result with Ultherapy's obvious lifting result only after one session. HIFU helps us attract more and more customers, love it!



Norman Bethune

The results are simply incredible specially in the area around my cheeks and mouth where I can see a constant smoothness and disappearance of fine and expression lines. Also the wrinkles around my eyes and forehead have been substantially minimised. You can see results in 3 weeks.

Before & After



Technical Specification

Irradiation Type

High Intensity Focused Ultrasound

Frequency

4MHz, 7MHz, 2MHz, 3MHz, 4MHz, 7MHz, 10MHz

Max Energy

2.0J/m²

LCD Screen

15 Inch True Color Touch LCD ADSS System

Ultra gun handle

10MHz 1.5mm; 7MHz 3mm; 4MHz 4.5mm

Face handle

10MHz 1.5mm; 7MHz 3mm; 4MHz 4.5mm

Body handle

3MHz 7mm; 2MHz 13mm

Voltage

AC110-220V 50/60Hz